



Computer Aided Cancer Diagnosis for serous effusion cytology

H. Elie¹, M. Lécluse¹, O. Lézoray², A. Elmoataz² ¹ Cherbourg Public Hospital Centre, Cherbourg, France ² University of Caen and CNRS, Caen, France





Serous effusion cytology


Types smears Papanicolaou using international stain:


- Conventional
- Cytospin
- Megafunal

ThinPrep









SCN400

Partial scanning with a regular sampling (200 to 500 images per region of interest): >100.000 un-compressed images 1024x1024

268 images with WSI scanner Aperio


249 images with WSI scanner Leica

517 images


40000x30000




50 Normal cytologies with only inflammation



20 Sub-normal cytologies with clinical suspicion of tumor





50 Adenocarcinoma metastasis




20 Mesothelioma

System benefits










- Gallery of cells, cells formula
- Cell thumbnail and context
- Feature statistics
 - For a single cell
 - For cell populations
 - Between two cell databases
 - Histograms and plots
- Easy remote expertise with several pathologists (review only suspicious cells).

Conclusion

- Digital cytopathology is a precious aid for pathologists in:
 - > Early diagnosis of tumors
 - Quality control
- Our system can be easily adapted to other cytologies
 - > Work in progress for the gynecological field


Image Segmentation/Analysis




Up to 100 features are computed per extracted cellular object:

- Size
- Shape
- Color
- Optical density Texture





An architecture of cascaded neural networks is used to classify extracted cellular objects:





Performance Evaluation

- 150 slides analyzed with certified tumor type: 98.5% of these slides recognized by our system (the 1.5% false negatives were due to out of focus problems).
- Cases with a clinical notion of tumor: 2.5% of cases diagnosed as negative by pathologists (3 out of 125) were evaluated as abnormal by our system (and it was really the case). The pathologist is thus alerted. For example the system extracted these rare abnormal cells:











- The system is efficient for rare event detection: on a specific case, 373 abnormal cells out of a total of 107,102 extracted cells were found (0.3%).
- 1000 cells sorted per second / processing time around 3 seconds per image (1024x1024 or 1360x1024 pixels)

Our Partners



